

FerleTest narkotika urintest

Tester for:

| <u>Test</u> | <u>Substans</u> | <u>Grænseværdi (ng/ml)</u> |
|-------------|---|----------------------------|
| THC | 11-nor- Δ^9 -THC-9-COOH | 200 |
| | 11-nor- Δ^8 -THC-9-COOH | 200 |
| | 11-hydroxy- Δ^9 -THC | 2,500 |
| | Δ^8 -Tetrahydrocannabinol | 7,500 |
| | Δ^9 -Tetrahydrocannabinol | 10,000 |
| | Cannabinol | 10,000 |
| | Cannabidiol | 100,000 |
| Tramadol | Tramadol | 200 |
| | N-desmethyl-tramadol | 500 |
| | O-desmethyl-tramadol | 20,000 |
| Buprenorfin | Buprenorphine-3- β -d-glucoronide | 10 |
| | Buprenorphine | 200 |
| Ketamin | Ketamine | 1,000 |
| | Norketamine | 500 |
| | Phencyclidine | 25,000 |
| | Tetrahydrozoline | 50,000 |
| | Chlopheniramine | 100,000 |
| | Dextromethorphan | 100,000 |
| | Lidocaine | 100,000 |
| | Promethazine | 100,000 |
| | d-Pseudoephedrine | 100,000 |
| Fentanyl | Fentanyl metabolites | 10 |

| <u>Test</u> | <u>Substans</u> | <u>Grænseværdi (ng/ml)</u> |
|---------------------------|---|----------------------------|
| Metamfetamin (Ecstasy) | (+)Methamphetamine | 1000 |
| | (±)3,4Methylenedioxyamphetamine (Ecstasy) | 1000 |
| | d-Amphetamine | > 100 µg/ml |
| | l-Amphetamine | > 100 µg/ml |
| | (±)3,4Methylenedioxyamphetamine | > 100 µg/ml |
| | Chloroquine | > 100 µg/ml |
| | (-)Ephedrine | > 100 µg/ml |
| | β-Phenylethylamine | > 100 µg/ml |
| | Procaine | > 100 µg/ml |
| | d-Pseudoephedrine | > 100 µg/ml |
| | Randitidin | > 100 µg/ml |
| Benzodiazepiner | Nitrazepam | 100 |
| | Chloridiazepoxide HCl | 300 |
| | Clobazam | 300 |
| | Desmethyldiazepam | 300 |
| | Oxazepam | 300 |
| | Temazepam | 300 |
| | Alprazolam | 1000 |
| | Bromazepam | 1000 |
| | Diazepam | 1000 |
| | Flunitrazepam | 1000 |
| | Lorazepam | 1000 |
| Clonazepam | 2000 | |
| Flurazepam | 100 | |
| Metadon | Methadone | 300 |
| | Methadol | 300 |

| <u>Test</u> | <u>Substans</u> | <u>Grænseværdi (ng/ml)</u> |
|---------------|-------------------------------------|----------------------------|
| Barbiturater | Alphenal | 100 |
| | Barbital | 150 |
| | Pentobarbital | 150 |
| | Phenobarbital | 150 |
| | Amobarbital | 300 |
| | Secobarbital | 300 |
| | Butalbital | 5,000 |
| Kokain | Benzoylecgonine | 300 |
| | Cocaine | 30,000 |
| Amfetamin | D-Amphetamine | 1,000 |
| | D/L-Amphetamine | 2,000 |
| | (±)3,4Methylenedioxyamphetamine | 2,500 |
| | l-Amphetamine | 30,000 |
| | (+)methamphetamine | > 100 µg/ml |
| | (±)3,4Methylenedioxymethamphetamine | > 100 µg/ml |
| Buprenorfin | Buprenorphine-3-β-d-glucoronide | 10 |
| | Buprenorphine | 200 |
| Phencyclidine | PCP | 25 |